



## Newsletter

You may remove your email address from our newsletter mailing list at anytime by using the opt out at the bottom of the page.


### Featured Product - The *newly redesigned* GCA-03(W) with electronic analog meter.

The analog meter is calibrated in mR/hr that is useful for approximating radiation levels. There are three pre-set ranges:

#### Ranges

- 0-1 mR/hr (Range switch set to 1X)
- 0-10 mR/Hr (Range switch set to 10X)
- 0-100 mR/hr (Range switch set to 100X)

The geiger counter has a digital output that provides a TTL logic pulse each time a radioactive particle is detected. Used to connected to various digital adapters and network software.



[More Information...](#)

### Geiger Counter Belt Clip

All of our assembled, in-case, geiger counters now have an optional add on belt clip. It is available on the following models:

[GCA-04/04W](#)

[GCA-03/03W](#)

[GCK-01A / GCK-02A](#)



[Become a fan of Kirlian Photography on Facebook!](#)

[Our new forum](#)

[Our website](#)

[HTML version](#)

[PDF version](#)

© 2011 Images SI, Inc.